

Introducing Your Next Monitor!

BM3Vet Next



Ph: (800) 211-0036



Fx: (888) 729-9010



www.apexx360.com

1. At a Glance

| | | |
|--|----------------------|-----------------|
| 7" High Resolution Widescreen (800 x 480) TFT | 2 ECG traces | 120 hours trend |
| LAN & HL7 | ST Analysis | PVC Count |
| Limited Arrhythmia Analysis | Large Parameter Mode | Trend Graph |
| Simplified User Interface | | |
| Still The Industrie's Longest Warranty at 4 Years! | | |



Ph: (800) 211-0036




Fx: (888) 729-9010




www.apexx360.com


2. Lots of Display Modes




Main Display Mode



ST Analysis Display Mode




Large Parameter Display Mode



Trend Display Mode



3. Display adjusts if functions are turned off



SpO2 function is turned off

