

apexx AE-1000 Portable X-Ray

The apexx AE-1000 Portable X-Ray is perfect for small animal or equine in or out of the clinic.

This portable x-ray's functionality and durable design provides years of productive use.



AC Input	VAC - 110VAC, 50/60 Hz
X-ray Tube	Toshiba
Focal Spot Size	1.8 mm
Total Inherent Filtration.....	0.7mm Al
Collimator	Dual knob, fully adjustable
Exposure Switch.....	Two stage
Exposure Cord	10 feet
Power Cord	16 feet
Size.....	9.5"W x 8.5"H x 10.5"L
Weight.....	37 lb (17 kg)

Output:

Electronic time selection; 23 steps available. Range: 0.03 – 5.0 seconds

mA / kVp Stages (5 total)	Time Range
10 mA @ 100 kVp	0.03 – 5.0 seconds
15 mA @ 90 kVp	0.03 – 5.0 seconds
20 mA @ 80 kVp	0.03 – 5.0 seconds
25 mA @ 70 kVp	0.03 – 5.0 seconds
30 mA @ 60 kVp	0.03 – 5.0 seconds

Every system includes: AE-1000 Portable X-Ray, power cord, exposure switch, carrying case, operators manual.

Contact Information:



apexx Veterinary Equipment
 14 Inverness Drive. East, Suite C-124
 Englewood, CO 80112

Ph: 800.211.0036
 Email: help@apexx360.com

www.apexx360.com

* Information and specifications subject to change without notice.