

apexx AE-1020HF High Frequency Portable X-Ray

The apexx AE-1020HF High Frequency Portable X-Ray is an efficient high frequency portable x-ray that is light weight, durable and powerful.

The small x-ray tube focal spot produces an excellent image with fine detail and superb resolution.



AC Input	VAC - 110VAC, 50/60 Hz
X-ray Tube	Toshiba
Focal Spot Size	1.2 mm
Anode Heat Storage	20 kHU
Total Filtration	2.6 mm Al
Collimator	Dual knob, fully adjustable
Exposure Switch.....	Two stage
Exposure Cord	9 feet
Power Cord	15 feet
Size.....	9"W x 7"H x 12"L
Weight.....	24.2 lb (11 kg)

Output:

kVp Range (2kV steps)	mAs Range
40 – 48 kVp	0.4 – 50 mAs
50 – 58 kVp	0.4 – 40 mAs
60 – 68 kVp	0.4 – 40 mAs
70 – 78 kVp	0.4 – 40 mAs
80 – 88 kVp	0.4 – 40 mAs
90 – 100 kVp	0.4 – 40 mAs

Every system includes: AE-1020HF High Frequency Portable X-Ray, power cord, exposure switch, carrying case, operators manual.

Contact Information:



apexx Veterinary Equipment
14 Inverness Drive. East, Suite C-124
Englewood, CO 80112

Ph: 800.211.0036
Email: help@apexx360.com

www.apexx360.com

* Information and specifications subject to change without notice.