

Sentinel-V 3020

Multi-Parameter Vital Signs Monitor

Patient Monitoring

The Sentinel-V 3020 provides more features and flexibility than Cardell/Midmark, Surgivet and Scil multi-parameter monitors, at a lower price.

The Ultrabright 8.4" display is viewable from any angle, providing clear view of waveforms and full monitoring parameters. Lightweight and portable, it's designed to monitor all sizes and types of animals, exotics to equine.

The 3020 multi-parameter monitor utilizes proven and dependable technology.

Both Nellcor OxiMax SpO2 technology and OxyCRG technology are standard and powerful. OxyCRG displays the interactive relation between HR, Resp and Oxygen, convenient for observation of clinical changes in smaller patients.



Product Highlights

- Ultrabright 8.4" display
- ECG / SpO2 / NIBP / Resp / Temp / Printer (EtCO2 optional)
- Nellcor OxiMax SpO2 technology
- 120-hour graphic and tabular trends of all parameters
- Built-in recorder and rechargeable battery

Sentinel-V 3020

Multi-Parameter Vital Signs Monitor

Patient monitoring is a must in every veterinary clinic. The Sentinel-V 3020 is a multi-parameter vital signs monitor that measures all standard parameters, viewable on an Ultrabright 8.4" display.



Patient Monitoring

Specifications

- Ultrabright 8.4" display
- ECG / SpO2 / NIBP / Resp / Temp / Printer (EtCO2 optional)
- Nellcor Oximax SpO2
- Li-ion battery provides 4 hours of operating time
- USB port for data transfer
- Portable and Lightweight
- BPM 20-320
- Up to 9 waveforms on one screen
- 120-hour graphic and tabular trends of all parameters
- HRV, Arrhythmia analysis and S-T segment analysis
- Resistant to high-frequency interference
- Built-in recorder and rechargeable battery

Standard Accessories:

- 3-lead ECG Lead Cable and Clips
- Y-Clip Lingual Sensor
- Assortment of five Cuffs 3cm - 17cm
- Paper

Optional Accessories:

- Mainstream EtCO2 Sensor
- Temperature Probe
- Esophageal Probes (3 sizes)
- 5-Lead ECG Cable and Clips
- IV Pole Clamp
- Mobile Stand

Leading Edge provides Veterinary professionals high-performance, competitively priced, reliable Veterinary medical equipment, backed by world-class support and service. Through our Green Initiative, we've made a genuine commitment to support long term environmental sustainability and we're committed to protecting and preserving the environment. As a single source Veterinary medical equipment company, Leading Edge products are available only through a limited and select group of Distributors, who share in our philosophies and standards of excellence, commitment and service.



Menu Knob

Push the Menu knob and operational menus will appear on the 8.4" display. Menu functions hi-light as you turn the knob. Simply select the desired function by pressing the knob.



Connection Ports

Color-coded, labeled ports make cable, probe and sensor connections simple.



Printer

Built-in recorder is standard.

Authorized Distributor

