

8 I G * \$ J 9 H I ` f U g c i b X

TRANSDUCER OPTIONS



Micro-convex Array Transducer*

Applications: Abdomen - small ; Obstetric - small ; Cardiac - small & large
Scanning Angle: 128 degrees
Scanning Depth: 157mm
(5.5/6.5/7.5/H4.5/H4.7 MHz)
Product ID: C611-2



Convex Array Transducer

Applications: Large animal: Abdomen, Obstetric
Scanning Angle: 60 degrees
Scanning Depth: 324mm
(2.5/3.5/4.5/H2.5/H2.7 MHz)
Product ID: C361-2



Endorectal Transducer

Applications: Large animal : Obstetric, tendon, back fat
Scanning Length: 60mm
Scanning Depth: 176mm
(4.5/5.5/6.5/H4.0/H4.2 MHz)
Product ID: V563-2



Linear Array Transducer

Applications: Small animal - Abdomen, Muscu-loskeletal, tendon, small parts, peripheral vascular.
Scanning Length: 40mm
Scanning Depth: 117mm
(6.5/7.5/8.5/H4.5/H4.7 MHz)
Product ID: L743-2



Linear Array Transducer

Applications: Small animal - Abdomen, Muscu-loskeletal, tendon, small parts, peripheral vascular.
Scanning Length: 60mm
Scanning Depth: 117mm
(6.5/7.5/8.5/H4.5/H4.7 MHz)
Product ID: L761-2



Mobile Trolley

Description: 18"W x 28"D x 36"H, 4" casters with brakes, accessory drawer and printer holder.
Product ID: DUS -Trolley

Imaging



DUS 60 VET Ultrasound includes :

Portable Ultrasound system for cardiology, peripheral vessels, abdomen, urology, pregnancy, small parts and biopsy assessment. Pseudo color PW Doppler, 12.1" LCD monitor, 256-frame cine loop memory, 504MB storage, 256 Gray scale, two USB ports, two transducer parts, Li-ion battery, gel and users manual. *INCLUDES: C611-2 Micro-convex Probe(5.5/6.5/7.5/H4.5/H4.7)

Warranty: Systems: 24 months; Transducers, Mobile Trolley, Footswitch; Sony Printer: 12 months; Storage/Transport Case and DICOM 3.0 software are warranted for out-of-box manufacturer defect.

Ultrasound Gel, Video Printer-B&W, paper, Storage/Transport Case, Footswitch, Handswitch, Needle Guide Bracket and DICOM 3.0 Software are also available.