

Sentinel-V 3025i

Multi-Parameter Vital Signs Monitor

Patient Monitoring

The Sentinel-V 3025i provides more features and flexibility than Cardell/Midmark and Surgivet multi-parameter monitors, at a lower price.

The 8.4" Touch Screen display is viewable from any angle, providing clear view of waveforms and full monitoring parameters.

The 3025i multi-parameter monitor features proven and dependable technology.

It now comes with patient data transfer capabilities, Nellcor OxiMax SpO₂, and first in industry OxyCRG and Adap-DSP™, providing advanced diagnostic information for all sizes and types of animals. These features make the best vital signs monitor even better.



Product Highlights

- 8.4" Touch screen display
- ECG / SpO₂ / NIBP / Resp / Printer / Temp & EtCO₂ optional
- Adap-DSP™ provides precise NIBP measurement
- Nellcor OxiMax SpO₂ technology
- Display freeze and ECG rewind
- Patient data transfer

Sentinel-V 3025i

Multi-Parameter Vital Signs Monitor

Patient monitoring is a must in every veterinary clinic. The Sentinel-V 3025i is a multi-parameter vital signs monitor that measures all standard parameters viewable on an 8.4" Touch Screen display.



Patient Monitoring

Specifications

- 8.4" Touch Screen display
- ECG / SpO2 / NIBP / Resp / Printer (Temp + EtCO2 sensors optional)
- Mainstream EtCO2 capable (Plug 'n Play option)
- Adap-DSP™ adaptive filter algorithm provides precise NIBP measurement
- Nellcor OxiMax SpO2
- Nellcor SatSeconds
- Patient data transfer
- Display Freeze and ECG Rewind
- Patient Data Capability: name, sex, weight, type
- Li-ion battery provides 4 hours of operating time
- Portable and Lightweight
- BPM 20-320
- Up to 7 waveforms on one screen
- 96-hour graphic and tabular trends of all parameters
- HRV, Arrhythmia analysis and S-T segment analysis
- Resistant to high-frequency interference
- Built-in printer and rechargeable battery

Standard Accessories:

- 3-lead ECG Lead Cable and Clips
- Y-Clip Lingual Sensor
- Assortment of five Cuffs 3cm - 17cm
- Paper

Optional Accessories:

- Mainstream EtCO2 Sensor
- Temperature Probe
- Esophageal Probes (3 sizes)
- 5-Lead ECG Cable and Clips
- IV Pole Clamp
- Mobile Stand

Leading Edge provides Veterinary professionals high-performance, competitively priced, reliable Veterinary medical equipment, backed by world-class support and service. Through our Green Initiative, we've made a genuine commitment to support long term environmental sustainability and we're committed to protecting and preserving the environment. As a single source Veterinary medical equipment company, Leading Edge products are available only through a limited and select group of Distributors, who share in our philosophies and standards of excellence, commitment and service.



Touch Screen Display

The 8.4" Touch Screen Display makes navigating the menu easier than ever.



Adap-DSP™

Combines quicker deflation times with adaptive signal processing to deliver more precise NIBP measurements



Nellcor SatSeconds

Reduces "nuisance" alarms for more reliable notifications

Authorized Distributor