

VE-H100B

Pulse Oximeter

Take a second look at the VE-H100B

The VE-H100B handheld Pulse Oximeter rivals the monitoring "leaders" for breadth of features, durability and value. The VE-H100B is specifically designed for veterinary use featuring veterinary specific software with enhanced sensitivity.

Advanced signal processing technology is ideal for patient motion and low perfusion cases and can be used for spot checks or continual monitoring during surgical and dental procedures.

Accurate Patient Monitoring in the Palm of Your Hand

The VE-H100B Pulse Oximeter is lightweight and portable and has a bright LCD display which shows SpO₂, Pulse Rate and waveforms. It offers up to 300 hours of data storage which can be transferred to a PC for storage, management, review, and printing.

Patient Monitoring



Product Highlights

- Handheld - easy to use
- Pulse rate range 30-350 bpm accomodates equine to exotics
- Sampling flexibility with FIVE available sensors
- Plethysmographic waveform/data transfer
- Long battery life - up to 48 hours on 4 "AA"




VE-H100B

Pulse Oximeter

Accurately Monitor Your Patients

Assessing oxygen saturation is essential in patient monitoring. The VE-H100B Pulse Oximeter is an accurate, affordable and easy-to-use pulse oximeter.



Display Screen

Large numeric and Waveform display with Real Time Clock, Perfusion Index (PI) Bar, Trend Graph, and SpO2 Pulse Review Table.



Easy-to-use Menu

Simple navigation buttons are conveniently located for ease of operation. Backlight switch and automatic power-off function saves power.



Available Accessories

Boot covers, Y-Clip, C-Clip, Rectal SpO2 Sensors, and Data Management Kit.

Patient Monitoring

Specifications

- LCD Display
- Real Time Clock (RTC) Display
- Backlight switch and automatic power-off function
- SpO2 and Pulse Rate Measurement
- Display number and waveform of SpO2 simultaneously
- Long battery life - up to 48 hours on 4 "AA"
- Powerful data storage capacity (up to 300 hours)
- Data transferable to PC for storage, review and printing
- Trend graph and table of SpO2 and Pulse review on the screen
- Audible and visual alarms
- Pitch tone

Standard Accessories:

- Y-Clip Lingual Sensor
- Small and Large Clips
- 4 x AA batteries
- Zippered Carrying Case

Optional Accessories:

- C-Clip SpO2 Sensor
- Rectal SpO2 Sensor
- Flexible Wrap SpO2 Sensor
- Data Management Kit
- Protective Boot Covers (Blue, Green, or Grey)

Leading Edge provides Veterinary professionals high-performance, competitively priced, reliable Veterinary medical equipment, backed by world-class support and service. Through our Green Initiative, we've made a genuine commitment to support long term environmental sustainability and we're committed to protecting and preserving the environment. As a single source Veterinary medical equipment company, Leading Edge products are available only through a limited and select group of Distributors, who share in our philosophies and standards of excellence, commitment and service.

Authorized Distributor



VER.072622