

Veterinary DR System




IWV System X

IWV Flat Panel Detector

- 17" x 17" Tethered or 14" x 17" Wireless
- Cesium CSi
- Pixels: 3072x3072 / Pixel Pitch: 140um
- Also available to upgrade most film tables





IWV System X

- 21" wide "Heads Up" Touchscreen display
- PC Workstation, Acquisition Software and Touchscreen Monitor
- Anatomical Programming (APR)
 - + Touch and Image
- 4-way float top table
- Hands-free collimator light control

High Frequency Generator

- 32kW / 400mA 125kVp
- Space-saving, inside table base
- Optional 110V power

Single Operator Heads-Up" Control Interface

- Four-step operation without stepping away from table or patient:
 - + Select APR (Touchscreen)
 - + Hands-free collimation
 - + Position patient (4-way float)
 - + Expose—Preview—Save
- Cuts personnel cost in half
- Image available in seconds = Immediate retakes
- Simple. Efficient. Convenient.
- Easier=More radiographs
- More revenue
- Better medicine