

# Multi-Purpose Table

## Lift & Transport Table

### Truly Multi-Purpose

Veterinarians utilize the Leading Edge Multi-Purpose Table as a treatment table, a vehicle to lift and transport patients in and out of the hospital, perform exams and surgery, and safely lift and transport heavy inventory.

### Reliable, Smooth & Quiet Hydraulic Lift

Easy to raise and lower, never needs plugged-in and always ready to work. Modular design for affordable parts replacement.

### What Your Colleagues Are Saying

*"We just received our 5th Leading Edge Multi-Purpose Table. We LOVE THESE! I am a customer for life. Don't know why more practices don't use these. I highly recommend buying one or a few."*

– C. McMillen, DVM

### Patient Management




### Product Highlights

- Reliable, smooth and quiet hydraulic lift
- Easily lifts and transports 300+ lbs.
- Lift range: 11" - 36"
- Balloon-tread casters roll smoothly
- Trailing-edge brakes are effortless to engage
- Steering handle allows safe maneuverability

# Multi-Purpose Table

## Lift & Transport Table



### Maintenance-free Hydraulic Lift

Never needs plugged-in.  
Always ready.

### Engineered to Provide Years of Service

Features include a large 22" x 48" stainless steel top with anti-drip edge. Its four corner bumpers protect your walls, doors and door frames.



### Balloon-tread Dual-Locking Casters

Roll smoothly and easily in and out of the clinic. Trailing-edge brakes are easy to engage and disengage, won't corrode or damage shoes or pants.



### User-Selectable Casters (Optional)

Provides additional control and safety; makes the Multi-Purpose Table easier and safer to move.



### Qwik-Set Tie Downs (Optional)

Technician's favorite option allows quick and easy patient restraint and release.

## Patient Management

### Specifications

- Lift range: 11" - 36"
- Trailing-edge brake tabs are easy to find, effortless to engage and disengage, won't rust or seize and won't tear up shoes or pant legs
- 22" x 48" Stainless steel top with anti-drop edge, drain hole and hook for pail
- Four corner bumpers protect your walls, doors and door frames
- Four tie-downs, two on each side
- Steering handle saves your back and provides easy and safe maneuverability without having to lean down to push or pull the table

#### Standard Accessories:

- 5" Dual-locking casters
- Standard tie-downs (4)

#### Optional Accessories:

- 6" Balloon-tread dual-locking casters
- User-selectable casters
- Qwik-set tie-downs
- Restraint straps

Leading Edge provides Veterinary professionals high-performance, competitively priced, reliable Veterinary medical equipment, backed by world-class support and service. Through our Green Initiative, we've made a genuine commitment to support long term environmental sustainability and we're committed to protecting and preserving the environment. As a single source Veterinary medical equipment company, Leading Edge products are available only through a limited and select group of Distributors, who share in our philosophies and standards of excellence, commitment and service.

### Authorized Distributor



VER.101222