

## Two Great Dental Stations, But What's the Difference?

### Prestige

### Ultimate

Asepsis, Modern Design .....	Traditional Design
2 Year Warranty .....	1 Year Warranty
LED Lighted Piezo Scaler .....	Non Lighted Piezo Scaler
Replaceable Scaler Handpiece - "Disconnectable" Handpiece .....	Must Replace Scaler Tubing Assembly to Replace Handpiece
Large Internal Storage Area Under Lid .....	No Internal Storage
Side Located Recessed Operation Controls .....	Controls Located on Front Panel
Enhanced "Service - Ability", All Service Parts are Easily Accessed .....	Delivery System With Service Parts Located Under Top counter top
The Industry's ONLY Pneumatic lift, easily lift and adjust height of unit for easy storage (under cut away wet table) or Technician operational comfort .....	Manual Height Adjustment
Oiless Compressor Available .....	Oiless Compressor Available

