

SUPERA

ANESTHESIA
INNOVATIONS

Safety Pressure Relief Valve

A revolutionary advance in anesthesia safety!



Features:

- Pressure relief at 20-25cmH2O using ventilation button
- Less than 1cmH2O pressure at full open
- Less than 1cmH2O pressure at full closed
- Stainless steel enclosure
- Fits any industry standard anesthesia machine with 3/4-20 thread
- Proven life saving patented design
- Can be used with ventilators

Part Number: ACC311



SUPERA

ANESTHESIA
INNOVATIONS

(503) 723-5068 or
(877) 620-1500 toll-free

info@superavet.com

12330 S.E. Highway 212
Clackamas, OR 97015